



April 12, 2018

2018 NNAAP® Equipment List for Skills Evaluation

Greetings:

Pearson VUE has revised the Equipment List for the NNAAP Skills Evaluation scheduled for administration beginning in July 2018.

The changes include the addition of XXL clothing tops for "Dresses Client with Affected (Weak) Right Arm," requirement that the PPE isolation gown secure in the back or the side (cannot secure in the front), and a small increase in the number of washcloths and towels required per candidate testing. Although the equipment list changes are not required until the administration of the revised Skills Evaluation scheduled for July, Pearson VUE is providing this preview to help with your planning.

If you have questions please contact <u>pearsonvuenurseaide@pearson.com</u> or call (877-538-3415).

Sincerely,

Jennifer Mankowski, BSN, RN

Jannifer Mankowski

Manager Nurse Aide Training and Auditing

Attachment: (1)

■ 2018 NNAAP® Skills Equipment List

EQUIPMENT LIST: NNAAP SKILLS EVALUATION AUDIT *JULY 2018



START-UP FURNIT				
ITEM	YES	NO	RECOMMENDED AMOUNTS	YES YES NO
sed, adjustable bed with non telescoping side alls (must be working). Head of the bed must go up/down and wheels lock/unlock			1 per NAE testing	
Chair			1 per NAE testing	
Clock (preferred) or Wrist Watch with Second Hand			1 per NAE testing unless wall clock can be clearly seen from all testing areas	· C
ax machine- Dedicated (accessible to skills lab)			1 per testing site	
inen Receptacle (designated for used/soiled inen)			1 per NAE testing	Or Similar to be Marked as "Dirty Linen"
Mannequin* (must be a full fernale mannequin able to accept catheter)*[CA: Although no catheter skill berformed, mannequin with female genitalia is required or perineal care,]*[WI only: partlal female mannequin that accepts catheter is acceptable]			1 per every 2 NAE testing	
Privacy curtain*, Screen or Door if private room *North Carolina- Must be a 360° Privacy curtain]			1 per NAE testing	Or Moveable Screens or Door
Scale, calibrated (bathroom/standing) NO electronic scale, NO digital scale] North Carolina- Must be a Standing Scale]			1 for usage of all NAEs	
Signaling device, movable (must be on cord and attached to bed- can be non-functional)			1 per skills bed	Found by or attached to hospital bed
Sink with running water in room (hot /cold)			minimum 1 sink	
Table, bedside			1 per NAE testing	
Table, over bed			1 per NAE testing	
Toilet/Bedside Commode clearly labeled toilet			1 for usage of all NAEs	F
Wastebasket with liner			1 by each bed used; 1 by each sink used	
Wheelchair with footrests			1 for usage of all NAEs	476

PREVIEW - EQUIPMENT LIST: NNAAP SKILLS EVALUATION AUDIT *JULY 2018

TEM YES NO RECOMMENDED AMOUNTS Salth towels 1 3 per Candidate Scheduled 1 2 per NAE testing 2 per NAE testing 3 per Candidate scheduled 1 2 per NAE testing 3 per Candidate scheduled 1 2 per NAE testing 3 per Candidate scheduled 2 per NAE testing 3 per Candidate testing 4 per Candidate testing 5 per NAE testing 5 per NAE testing 5 per NAE testing 6 per NAE testing 6 per NAE testing 7 per NAE testing 7 per NAE testing 8 per NAE testing 9 per NAE testing 1 per skills bed 1 per skills bed 3 per NAE testing 4 per Candidate scheduled 5 per NAE testing 6 per NAE testing 1 per skills bed 3 per NAE testing 1 per skills bed 3 per NAE testing 4 per Candidate scheduled 5 per NAE testing 6 per NAE testing 6 per NAE testing 1 per skills bed 3 per NAE testing 4 per Candidate testing 5 per NAE testing 6 per NAE testing 1 per skills bed 3 per NAE testing 4 per Candidate testing 5 per NAE testing 6 per NAE testing 7 per NAE testing 8 per NAE testing 9 per NAE testing 1 per NAE tes	LINEN	CLOT!	IING	1997年1997年1997年1997年	
Service (patient)	ITEM			RECOMMENDED AMOUNTS	
Soling protector (bits or towed) 2 per NAE testing 2 per NAE				PRINCE TO CONTROL OF BURNESS OF A CONTROL OF THE PROPERTY OF T	
Zeper NAE testing					
Joseph the Front, no more than our (e) (buttons/nraps -no hospital gowns) Johns bottom sheets (fitted or flat) Jeps NAE testing Geach per NAE testing Geach per NAE testing Jeps NAE testin					
Courte (Stational/Anapa-no-hospital gowns Sizewas (patient) 2 per candidate testing 3 per ANE testing 6 each per NAE testing 1 per ANE testing 1 per ANE testing 1 per ANE testing 2 per ANE testing 3 per ANE testing 4 per candidate scheduled 5 per ANE testing 4 per candidate scheduled 5 per ANE testing 6 per ANE testing 7 per ANE testing 7 per ANE testing 8 per ANE testing 9 per ANE testing 1 per ANE testing 2 per ANE testing 2 per ANE testing 3 per ANE testing 1 per ANE testing 1 per ANE testing 2 per ANE testing 3 per ANE testing 4 per and testing 5 per ANE testing 5 per ANE testing 1 per ANE testing 1 per ANE testing 2 per ANE testing 3 per ANE testing 4 per and testing 5 per ANE testing 6 per ANE testing 7 per ANE testing 7 per ANE testing 7 per ANE testing 8 per ANE testing 9 per ANE testing 9 per ANE testing 1				2 per NAE testing	
Jones (patient) 2 per rAME testing 3 per NAE testing 6 each per NAE testing 6 each per NAE testing 7 per NAE testing 8 per skille bed 9 per NAE testing 10 per skille bed 9 per NAE testing 10 per skille bed 9 per NAE testing 10 per skille bed 10 pe	[[] [] [] [] [] [] [] [] [] [
inersis bottom sheets (fitted or flat) 2 per NAE testing	four (4)buttons/snaps -no nospital gowns				
iners bottom sheets (fitted or flat) 2 per NAE testing	Gowns (patient)			2 per candidate testing	0 0
Jack/Linen Protector (may use towel or fravwheet as pad) Interest pillows as pad) Japer skills bed				102	7 00 1
Jack/Linen Protector (may use towel or fravwheet as pad) Interest pillows as pad) Japer skills bed					11). A
Jack/Linen Protector (may use towel or fravwheet as pad) Interest pillows as pad) Japer skills bed					
Pad/Unen Protector (may use towel or fravahete as pad) Illiand I	Linens: bottom sheets (fitted or flat)			2 per NAE testing	
Pad/Unen Protector (may use towel or fravahete as pad) Illiand I	Linens: nillowcase, top sheet			6 each per NAE testing	
drawshed tas pad) Jeer skills bed	Emons, pinotrouse, top emost				
drawshed tas pad) Jeer skills bed					
drawshed tas pad) Jeer skills bed					
Subportive devices (pillows, blanket rolls, verdiges) Washoloth SUPPLES Washoloth SUPPLES TIEM YES NO RECOMMENDED AMOUNTS Basin, smeals Bedpan* (standard) **[PA not required] Electronic Blood Pressure Monitors with Conference of the Conference	Pad/Linen Protector (may use towel or			2 per NAE testing	
Supportive devices (pillows, blanket rolls, wedges) Masholoth SUPPLIS ITEM YES NO RECOMMENDED AMOUNTS Basin, bath SUPPLIS ITEM YES NO RECOMMENDED AMOUNTS 4 per NAE testing Basin, bath 2 per NAE testing Basin, emeris Bedpan* (standard) *IPA not required] 2 per NAE testing 2 per NAE testing 3 per NAE testing 3 per NAE testing 4 per NAE testing 5 per NAE testing 6 per NAE testing 6 per NAE testing 7 per NAE testing 8 per NAE testing 9 per NAE testing 1 per NAE testing 2 per NAE testing 1 per NAE testing 1 per NAE testing 2 per NAE testing 1 per NAE testing 2 per NAE testing 3 per NAE testing 4 per NAE testing 5 per NAE testing 6 per NAE testing 6 per NAE testing 7 per NAE testing 8 per NAE testing 9 per NAE testing 9 per NAE testing 1 per NAE testing 2 per NAE testing 3 per NAE testing 4 per NAE testing 5 per NAE testing 6 per NAE testing 9 per NAE testing 1 per NAE testing 2 per NAE testing 3 per NAE testing 4 per NAE testing 5 per NAE testing 6 per NAE testing 7 per NAE testing 9 per NAE testing 1 per NAE testing 1 per NAE testing 1 per NAE testing 1 per NAE testing 2 per NAE testing 3 per NAE testing 4 per NAE testing 5 per NAE testing 6 per NAE testing 7 per NAE testing 9 per NAE testing 1 per NAE test	drawsheet as pad)			I I	
washoloth SUPPLIES ITEM YES NO RECOMMENDED AMOUNTS Basin, amesis Bedin, emesis Bedin, emesis	Pillows				
Supplies				3 per NAE testing	
TITEM YES NO RECOMMENDED AMOUNTS	wedges)			A per candidate scheduled	
Assin, bath 4 per NAE testing 2 per NAE testing 2 per NAE testing 3 per NAE testing 4 per NAE testing 5 per NAE testing 5 per NAE testing 5 per NAE testing 6 per NAE testing 6 per NAE testing 7 per NAE testing 8 per NAE testing 8 per NAE testing 8 per NAE testing 8 per NAE testing 9 per NAE testing 1 per NAE testing 2 per NAE testing 1 per NAE testing 2 per NAE testing 1 per NAE testing 2 per NAE testing 2 per NAE testing 3 per NAE testing 4 per NAE testing 4 per NAE testing 5 per NAE testing 5 per NAE testing 5 per NAE testing 6 per NAE testing 7 per NAE testing 7 per NAE testing 8 per NAE testing 9 per NAE testing 1 per NAE testing 1 per NAE testing 9 per NAE testing 1 per NAE testing		IPPLIES		4 per candidate scheduled	THE THE PERSON WAS ASSESSED.
Basin, bath 4 per NAE testing 2 per NAE testing 2 per NAE testing 2 per NAE testing 3 per NAE testing 4 monitors with Cuffs (Co, DoC, SC, NH Olity). Please see BP Monitor Specification Document (not brand specific) 4 monitor Specifications Document (not brand specific) 5 monitor Specifications Document (not brand specific) 6 monitor Specifications Document (not brand specific) 7 monitor Specifications Document (not brand specific) 8 stating and specific) 8 stating 11 No. Off or trage of all NAEs. 1 double a praphic stating stationscope per NAE testing 9 per NAE testing 1 per NAE testing 1 per NAE testing 1 per NAE testing 1 per NAE testing Denture Brush (toothbrush labeled for denture skill) 1 per NAE testing 2 per NAE testing 1 per NAE testing 1 per NAE testing Dentures 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large, must not use older type with loose flap in toe area of the properties of the pro				RECOMMENDED AMOUNTS	
Basin, emeds 2 per NAE testing 2 per NAE testing 2 per NAE testing 2 per NAE testing 3 per NAE testing 4 per NAE testing 5 per NAE testing 6 per arm molitors with Country Please see 8P molitors with monitors, Monitor specifical per NAE testing at Stathocope*, dual expice of mobile stand molitors portable, val mounted, or on mobile stand at Stathocope*, dual expice of manual 8P (AL, AX, CA, GA, Gaum, Mo, Miss, IX, CA), For an additional per NAE testing 1 XL curif for usage of all NAEs. 1 dual expices stathocope by a curification of the stand of the st					
Bedpan* (standard) *[PA not required] 2 per NAE testing Minimum 2 machines; Upper arm monitors/cuffs only - no wrist monitors, fulf on his manual monitors, fulf only - no wrist monitors, fulf on his manual monitor					
Bedpan* (standard) *[PA not required] 2 per NAE testing Minimum 2 machines; Upper arm monitors/cuffs only - no wrist monitors, fulf on his manual monitors, fulf only - no wrist monitors, fulf on his manual monitor					
Electronic Blood Pressure Monitors with Curfs**(Co, DoC, SC, NH ONLY). Please see BP Monitor Specific) Minimum 2 machines; Upper arm monitors/curfs only - no wrist monitors; Monitor Specificians Document (not brand specific) Manual Blood Pressure Curfs**(standard or extra large curfs) and Set Brick Stocope*, dual earpice of or Manual Blood Pressure Curfs**(standard or extra large curfs) and BP**(AL, AK, Ca, Ga, Guam, MD, Miss, MC, ND, PA, RI, TX, YM, ANM, UY, ONLY) Set Brick Standard curf per NAE testing 1 XL curf for usage of all NAEs, 1 dual earpices stethoscope per NAE testing 1 per Mannequin 1 per Mannequin 1 per MAE testing Several drops per NAE testing Denture Brush (toothbrush labeled for denture skill) 2 per NAE testing Denture cup w/lid 1 per NAE testing Dentures 1 set or 2 uppers or 2 lowers for usage of all NABs Coloring to simulate urine (for urinary output skill) Dentures 1 set or 2 uppers or 2 lowers for usage of all NABs Open Toe or Closed toe OK Knee-high elastic stockings (clean, large or extra large, must not use older type with loose flap in toe area) 1 set or 2 uppers or 2 lowers for usage of all NABs Open Toe or Closed toe OK Measuring container** (graduate cylinder or hal- at least 250 ml**s/cc*s units clearly visible)(no urinals)* Pen or pencil 1 of each available	Basin, emesis			2 per NAE testing	
Electronic Blood Pressure Monitors with Curfs**(Co, DoC, SC, NH ONLY). Please see BP Monitor Specific) Minimum 2 machines; Upper arm monitors/curfs only - no wrist monitors; Monitor Specificians Document (not brand specific) Manual Blood Pressure Curfs**(standard or extra large curfs) and Set Brick Stocope*, dual earpice of or Manual Blood Pressure Curfs**(standard or extra large curfs) and BP**(AL, AK, Ca, Ga, Guam, MD, Miss, MC, ND, PA, RI, TX, YM, ANM, UY, ONLY) Set Brick Standard curf per NAE testing 1 XL curf for usage of all NAEs, 1 dual earpices stethoscope per NAE testing 1 per Mannequin 1 per Mannequin 1 per MAE testing Several drops per NAE testing Denture Brush (toothbrush labeled for denture skill) 2 per NAE testing Denture cup w/lid 1 per NAE testing Dentures 1 set or 2 uppers or 2 lowers for usage of all NABs Coloring to simulate urine (for urinary output skill) Dentures 1 set or 2 uppers or 2 lowers for usage of all NABs Open Toe or Closed toe OK Knee-high elastic stockings (clean, large or extra large, must not use older type with loose flap in toe area) 1 set or 2 uppers or 2 lowers for usage of all NABs Open Toe or Closed toe OK Measuring container** (graduate cylinder or hal- at least 250 ml**s/cc*s units clearly visible)(no urinals)* Pen or pencil 1 of each available					
Electronic Blood Pressure Monitors with Curfs**(Co, DoC, SC, NH ONLY). Please see BP Monitor Specific) Minimum 2 machines; Upper arm monitors/curfs only - no wrist monitors; Monitor Specificians Document (not brand specific) Manual Blood Pressure Curfs**(standard or extra large curfs) and Set Brick Stocope*, dual earpice of or Manual Blood Pressure Curfs**(standard or extra large curfs) and BP**(AL, AK, Ca, Ga, Guam, MD, Miss, MC, ND, PA, RI, TX, YM, ANM, UY, ONLY) Set Brick Standard curf per NAE testing 1 XL curf for usage of all NAEs, 1 dual earpices stethoscope per NAE testing 1 per Mannequin 1 per Mannequin 1 per MAE testing Several drops per NAE testing Denture Brush (toothbrush labeled for denture skill) 2 per NAE testing Denture cup w/lid 1 per NAE testing Dentures 1 set or 2 uppers or 2 lowers for usage of all NABs Coloring to simulate urine (for urinary output skill) Dentures 1 set or 2 uppers or 2 lowers for usage of all NABs Open Toe or Closed toe OK Knee-high elastic stockings (clean, large or extra large, must not use older type with loose flap in toe area) 1 set or 2 uppers or 2 lowers for usage of all NABs Open Toe or Closed toe OK Measuring container** (graduate cylinder or hal- at least 250 ml**s/cc*s units clearly visible)(no urinals)* Pen or pencil 1 of each available	Podpan* (standard) *[PA not required]			2 per NAF testing	
Cuffs*(CO, DOC, SC, NH ONLY). Please see B Monitor Specifications Document (not brand specific) Manual Blood Pressure Cuffs*(standard or extra large cuffs) and Stethoscope*, dual earpice *for monitors, Monitor portable, wall mounted, or on mobile stand I standard cuff per NAE testing I st. cuff for usage of all NAEs. I dual earpice catheter skill) * (CA, MD: not required) I per Mannequin I per NAE testing Coloring to simulate urine (for urinary output skill) Denture Brush (toothbrush labeled for denture skill) Denture cup w/lid I per NAE testing Several drops per NAE testing Denture cup w/lid I per NAE testing Open Toe or Closed toe OK Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) Massuring container* (graduate cylinder or hat- at least 250 ml 's/cc''s units clearly visible)(no urinals)* Pen or pencil I of each available Testing only no wrist monitors, Monitor por passes and provided and provided and provided in the standard of the standard o	Bedpair (Standard) [FA not required]			L per tivit testing	
Cuffs*(CO, DOC, SC, NH ONLY). Please see B Monitor Specifications Document (not brand specific) Manual Blood Pressure Cuffs*(standard or extra large cuffs) and Stethoscope*, dual earpice *for monitors, Monitor portable, wall mounted, or on mobile stand I standard cuff per NAE testing I st. cuff for usage of all NAEs. I dual earpice catheter skill) * (CA, MD: not required) I per Mannequin I per NAE testing Coloring to simulate urine (for urinary output skill) Denture Brush (toothbrush labeled for denture skill) Denture cup w/lid I per NAE testing Several drops per NAE testing Denture cup w/lid I per NAE testing Open Toe or Closed toe OK Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) Massuring container* (graduate cylinder or hat- at least 250 ml 's/cc''s units clearly visible)(no urinals)* Pen or pencil I of each available Testing only no wrist monitors, Monitor por passes and provided and provided and provided in the standard of the standard o					
Monitor Specifications Document (not brand specific) I standard cutf per NAE sting I studing rouge of all NAEs. 1 dual earpiece of all NAEs. 1 dual earpiece stethoscope per NAE testing I per NAE testing Denture Brush (toothbrush labeled for denture skill) I per NAE testing Denture cup w/lid I per NAE testing I set or 2 uppers or 2 lowers for usage of all NABs Coloring to simulate urine (for urinary output skill) Dentures I set or 2 uppers or 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) I stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) Measuring container* (graduate cylinder or hat-at least 250 ml*s/cc*s units clearly visible)(no urinals)* Pen or pencil Monitor Specific value are please to make a standard cutf for usage of all NABs Adaptive for usage of all NABs I stocking per each Candidate that has Knee-High Skill (PA/WI) I per NAE testing Open Toe or Closed toe OK Cylinder or Hat OK	Electronic Blood Pressure Monitors with				
Istandard cuff per NAE Istandard or extra Istandard or extra Istandard or extra Istandard cuff per NAE Istandard or extra Istandard stethoscope*, dual earpice extra Istandard cuff per NAE Istandard cuff	Cuffs*(CO, DOC, SC, NH ONLY). Please see BP				
Manual Blood Pressure Cuffs* (standard or extra large cuffs) and Stethoscope*, dual earpiece *for Manual Blood Pressure Cuffs* (standard or extra large cuffs) and Stethoscope*, dual earpiece *for Manual Br* (AL, AK, CA, GA, Gum, MD, Miss, NC, ND, PA, RI, RV, TW, AM, MI, VI, ONLY) Catheter and Syringe* (use with mannequin for catheter skill) * 1 per NAE testing Denture Brush (toothbrush labeled for denture skill) * 1 per NAE testing Denture cup w/lid * 1 set or 2 uppers or 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) The stocking per each Candidate that has Knee-High skill (PA/WI) Measuring containe** (graduate cylinder or hat-at least 250 m1*s/cc's units clearly visible)(no urinals)** Pen or pencil * 1 of each available * 1 of each each candidate that has the each candidate	The state of the s				
large cuffs) and Stethoscope*, dual earpiece *for Manual 8* "(AL, AK, CA, GA, Guarn, MD, Miss, NC, ND, PA, RI, TK, YT, WA, NM, IV, ONLY) Catheter and Syringe* (use with mannequin for catheter skill) * [CA, MD: not required] Denture Brush (toothbrush labeled for denture skill) Denture Brush (toothbrush labeled for denture skill) Denture cup w/lid Denture cup w/lid Denture cup w/lid 1 per NAE testing 1 set or 2 uppers or 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) Measuring container* (graduate cylinder or hat- at least 250 ml *s/cc* s units clearly visible)(no urinals)* Pen or pencil Lesting 1 Xt. cuff for usage of all NABs (PA/MI) 1 per NAE testing 1 per NAE testing Open Toe or Closed toe OK Cylinder or Hat. OK	specific)			mobile stand	
large cuffs) and Stethoscope*, dual earpiece *for Manual 8* "(AL, AK, CA, GA, Guarn, MD, Miss, NC, ND, PA, RI, TK, YT, WA, NM, IV, ONLY) Catheter and Syringe* (use with mannequin for catheter skill) * [CA, MD: not required] Denture Brush (toothbrush labeled for denture skill) Denture Brush (toothbrush labeled for denture skill) Denture cup w/lid Denture cup w/lid Denture cup w/lid 1 per NAE testing 1 set or 2 uppers or 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) Measuring container* (graduate cylinder or hat- at least 250 ml *s/cc* s units clearly visible)(no urinals)* Pen or pencil Lesting 1 Xt. cuff for usage of all NABs (PA/MI) 1 per NAE testing 1 per NAE testing Open Toe or Closed toe OK Cylinder or Hat. OK	Manual Bland Brossura Cuffs#letandard or outra			1 standard cuff per NAE	
Manual 8P *[AL, AK, CA, GA, Guam, MD, Miss, NC, ND, PA, RI, TX, VT, WA, NMI, VI - ONLY] Catheter and Syringe* (use with mannequin for catheter skill) * [CA, MD: not required]					
Catheter and Syringe* (use with mannequin for catheter skill) 1 per MAE tsting Coloring to simulate urine (for urinary output skill) 2 penture cup w/lid 2 per NAE testing Denture cup w/lid 2 per NAE testing 1 per NAE testing 1 per NAE testing 2 per NAE testing 1 per NAE testing 2 per NAE testing 1 per NAE testing 1 per NAE testing 2 per NAE testing 3 per NAE testing 4 per NAE testing 5 per NAE testing Copen Toe or Closed toe OK Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) Measuring container* (graduate cylinder or hat- at least 250 ml *s/cc'*s units clearly visible)(no urinals)* Pen or pencil 1 per NAE testing 1 per NAE testing 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) 1 per NAE testing Open Toe or Closed toe OK Cylinder or Hat OK OK 1 per NAE testing 1 per NAE testing Or similar					
catheter skill) * [CA, MD: not required] 1 per MAE testing Coloring to simulate urine (for urinary output skill) 2 penture cup w/lid 2 per NAE testing 1 per NAE testing 1 per NAE testing 1 per NAE testing 2 penture cup w/lid 2 penture cup w/lid 3 per NAE testing 1 per NAE testing 1 per NAE testing 2 pentures 1 set or 2 uppers or 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) Measuring container* (graduate cylinder or hat-at least 250 ml *s/cc* s units clearly visible)(no urinals)* Pen or pencil 1 per NAE testing 1 per NAE testing 1 per NAE testing 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) 1 per NAE testing Open Toe or Closed toe OK Cylinder or Hat OK OK 1 per NAE testing	Control & Control Cont			stethoscope per NAE testing	LEON.
Denture Brush (toothbrush labeled for denture skill) 1 per NAE testing Several drops per NAE testing 1 per NAE testing 1 per NAE testing 1 per NAE testing 1 per NAE testing Open Toe or Closed toe OK Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe large; must not use older type with loose flap in toe large; must not use older type with loose flap in toe large.) Measuring container* (graduate cylinder or hat-at least 250 ml 's/cc 's units clearly visible)(no urinals)* Pen or pencil 1 per NAE testing Open Toe or Closed toe OK Cylinder or Hat OK OK 1 per NAE testing 1 per NAE testing Open Toe or Closed toe OK 1 per NAE testing Or similar					
Coloring to simulate urine (for urinary output skill) Denture cup w/lid 1 per NAE testing 1 set or 2 uppers or 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) Measuring container* (graduate cylinder or hat-at least 250 ml 's/cc's units clearly visible)(no urinals)* Pen or pencil 1 per NAE testing Open Toe or Closed toe OK Cylinder or Hat OK 1 per NAE testing 1 or similar	catheter skill) * [CA, MD: not required]			1 per Mannequin	
Coloring to simulate urine (for urinary output skill) Denture cup w/lid 1 per NAE testing 1 set or 2 uppers or 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) Measuring container* (graduate cylinder or hat-at least 250 ml 's/cc's units clearly visible)(no urinals)* Pen or pencil 1 per NAE testing Open Toe or Closed toe OK Cylinder or Hat OK 1 per NAE testing 1 or similar					Enter Company of the
Denture cup w/lid 1 per NAE testing 1 set or 2 uppers or 2 lowers for usage of all NABs Cylinder or 1 clean stocking servate large; must not use older type with loose flap in toe area) Measuring container* (graduate cylinder or hat- at least 250 ml 's/cc 's units clearly visible)(no urinals)* Pen or pencil 1 per NAE testing 1 set or 2 uppers or 2 lowers for usage of all NABs 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) 1 per NAE testing Open Toe or Closed toe OK Cylinder or Hat OK OK 1 of each available or similar	Denture Brush (toothbrush labeled for denture skill)			1 per NAE tsting	
Denture cup w/lid 1 per NAE testing 1 set or 2 uppers or 2 lowers for usage of all NABs Cylinder or 1 clean stocking servate large; must not use older type with loose flap in toe area) Measuring container* (graduate cylinder or hat- at least 250 ml 's/cc 's units clearly visible)(no urinals)* Pen or pencil 1 per NAE testing 1 set or 2 uppers or 2 lowers for usage of all NABs 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) 1 per NAE testing Open Toe or Closed toe OK Cylinder or Hat OK OK 1 of each available or similar					
Denture cup w/lid 1 per NAE testing 1 set or 2 uppers or 2 lowers for usage of all NABs Cylinder or 1 clean stocking servate large; must not use older type with loose flap in toe area) Measuring container* (graduate cylinder or hat- at least 250 ml 's/cc 's units clearly visible)(no urinals)* Pen or pencil 1 per NAE testing 1 set or 2 uppers or 2 lowers for usage of all NABs 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) 1 per NAE testing Open Toe or Closed toe OK Cylinder or Hat OK OK 1 of each available or similar					
Dentures 1 set or 2 uppers or 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) Measuring container* (graduate cylinder or hat- at least 250 ml 's/cc 's units clearly visible)(no urinals)* Pen or pencil 1 per NAE testing 1 set or 2 uppers or 2 lowers for usage of all NABs 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) 1 per NAE testing Open Toe or Closed toe OK Cylinder or Hat OK OK 1 of each available	Coloring to simulate urine (for urinary output skill)			Several drops per NAE testing	
Dentures 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) 1 set or 2 uppers or 2 lowers for usage of all NABs 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) Measuring container* (graduate cylinder or hat- at least 250 ml 's/cc 's units clearly visible)(no urinals)* Pen or pencil 1 set or 2 uppers or 2 lowers for usage of all NABs 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) To Hat OK 1 of each available or similar					
Dentures 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) 1 set or 2 uppers or 2 lowers for usage of all NABs 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) Measuring container* (graduate cylinder or hat- at least 250 ml 's/cc 's units clearly visible)(no urinals)* Pen or pencil 1 set or 2 uppers or 2 lowers for usage of all NABs 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) To Hat OK 1 of each available or similar			-		
Dentures 2 lowers for usage of all NABs Knee-high elastic stockings (clean, large or extra large; must not use older type with loose flap in toe area) 1 set or 2 uppers or 2 lowers for usage of all NABs 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) Measuring container* (graduate cylinder or hat- at least 250 ml 's/cc 's units clearly visible)(no urinals)* Pen or pencil 1 set or 2 uppers or 2 lowers for usage of all NABs 1 stocking per Mannequin or 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) To Hat OK 1 of each available or similar					
2 lowers for usage of all NABs Open Toe or Closed toe OK	Denture cup w/lid			1 per NAE testing	
2 lowers for usage of all NABs Open Toe or Closed toe OK					11572-15
2 lowers for usage of all NABs Open Toe or Closed toe OK					
Cylinder or Hat Open or Pencil Cylinder or similar	Dentures			1 set or 2 uppers or	
The search of the static stockings (clean, large or extra large; must not use older type with loose flap in toe area)	Delicales		1	2 lowers for usage of all NABs	Open Top or Closed top OK
large; must not use older type with loose flap in toe area) 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) Measuring container* (graduate cylinder or hat- at least 250 ml's/cc's units clearly visible)(no urinals)* Pen or pencil 1 clean stocking per each Candidate that has Knee-High Skill (PA/WI) 1 per NAE testing 1 of each available or similar	Vnog high elastic stockings (clean large or extra	-		1 stocking per Mannequin or	Open foe of closed toe on
area) Measuring container* (graduate cylinder or hat- at least 250 ml 's/cc 's units clearly visible)(no urinals)* Pen or pencil 1 of each available Cylinder or Hat OK OK 1 per NAE testing 1 of each available					
Measuring container* (graduate cylinder or hat- at least 250 ml 's/cc 's units clearly visible)(no urinals)* Pen or pencil 1 per NAE testing 1 of each available or similar		1			
Pen or pencil 1 of each available 1 of each available	200.00		1		
Pen or pencil 1 of each available or similar			1	1 per NAE testing	OK.
		-	-		
Transfer (gait) belt 1 per NAE testing	Pen or pencil			1 of each available	or similar
Transfer (gait) belt 1 per NAE testing					
	Transfer (goit) holt	1	1	1 per NAE testing	
	Transfer (gait) beit	1	1		

PREVIEW - EQUIPMENT LIST: NNAAP SKILLS EVALUATION AUDIT *JULY 2018

DISPOSA	BLE SU	PPLIES		
ITEM	YES	NO	RECOMMENDED AMOUNTS	
Alcohol wipe or Alcohol and cotton ball			2 per candidate scheduled	
Antimicrobial Spray/Wipes- clearly labeled			1 container per NAE testing	
Drinking cup (disposable)			8 per NAE testing	W.
Food (two varieties that can typically be fed in bite - size portions with a spoon) and beverage water			2 for each NAE testing	
Gloves, large, disposable non-latex			3 pr per candidate scheduled	
Gown,* Isolation (long sleeve that secure at neck and waist at back or side, cannot secure in front. Cloth or disposable)			2 per NAE testing	
Hand Sanitizer			1 container per NAE testing	
Hand Wipes (may use washcloth)			2 per candidate scheduled	
Lotion in pump container (hypoallergenic & unscented)			1 container per NAE testing	
Meal tray with "name" card			1 per Meal Tray	
Napkins (may use paper towels)			4-6 per NAE testing	
Paper plates			8 per NAE testing	
Paper towels			2 rolls for usage of all NAEs	
Plastic bags (for wastebasket)			2 per each wastebasket	
Soap, liquid in pump container (hypoallergenic &			1 container per NAE testing	
unscented/not rinseless/no foaming soaps)			12	
Spoons (disposable)			8 per NAE testing	
Toilet tissue or Wipes			1 roll per NAE testing or	
			2 wipes per candidate scheduled	
Toothbrush (individually wrapped)			4-6 per NAE testing	
Toothpaste (separate tubes labeled for mouth care			1 per NAE testing	
and for denture care)				